



**AZUR<sup>®</sup> CX**

Peripheral Coil System

OPTIMIZE  
OUTCOMES  
IN **MORE**  
**PROCEDURES**



**Available in 2 and 3 mm sizes**  
for small vessel embolization

**TERUMO**  
INTERVENTIONAL  
SYSTEMS

# AZUR<sup>®</sup> CX

Peripheral Coil System

## OPTIMIZE EMBOLIZATION

A unique balance of coil design and Hydrogel technology for a wide range of procedures


### Soft, flexible Hydrogel for efficiency and controlled delivery

- Superior volume and packing density<sup>1,2</sup>
- Sustainable, natural tissue proliferation—may reduce incidence of recanalization<sup>3,4</sup>
- Mechanical occlusion—less reliance on thrombus formation
- Soft feel with up to 30 minutes repositioning time



### Designed to form a solid core

Enables expansion between the gaps with hydrogel that will not be absorbed by the body<sup>3,4</sup>



### Ease of Deployment

- Use the AZUR<sup>®</sup> Detachment System for precise positioning and placement
- Minimize catheter manipulation with unique coil design



## PUSHING BOUNDARIES

Terumo Interventional Systems is **committed to your success** with innovative procedural solutions and ongoing support for your most challenging cases.

We are relentlessly seeking new ways to help you apply effective solutions and achieve **better outcomes for more patients.**



**TERUMO**  
INTERVENTIONAL  
SYSTEMS

# AZUR<sup>®</sup> CX

Peripheral Coil System

## OPTIMIZE EMBOLIZATION

A unique balance of coil design and Hydrogel technology for a wide range of procedures


### Soft, flexible Hydrogel for efficiency and controlled delivery

- Superior volume and packing density<sup>1,2</sup>
- Sustainable, natural tissue proliferation—may reduce incidence of recanalization<sup>3,4</sup>
- Mechanical occlusion—less reliance on thrombus formation
- Soft feel with up to 30 minutes repositioning time



### Designed to form a solid core

Enables expansion between the gaps with hydrogel that will not be absorbed by the body<sup>3,4</sup>



### Ease of Deployment

- Use the AZUR<sup>®</sup> Detachment System for precise positioning and placement
- Minimize catheter manipulation with unique coil design



## EXPAND YOUR OPTIONS

The availability of 2 and 3 mm sizes enables AZUR<sup>®</sup> CX to support more procedures, including:


### Small Vessel Embolization

- Gastrointestinal bleeds
- Pre Y-90 gastric embolization
- Prostate artery embolization
- Other pelvic vasculature embolizations

### Large Vessel Embolization

- Internal iliac procedures
- Pelvic congestion syndrome
- Varicocele embolization
- Peripheral aneurysms
- Arteriovenous malformations

Treat more patients with the potential to use fewer coils<sup>5</sup>



**Cross coverage**—  
Designed to fill vessel with no gap in center

**Solid core**—Supports neointima growth with Hydrogel technology

**Anchor**—Use for control in high-flow areas

**Soft and flexible**—  
Use directly out of the package

## TAKE CONTROL

Choose from a comprehensive line of peripheral coils from TERUMO\*

Embolization Goal	Aneurysms		Vessel occlusions	
	Microcatheter	Diagnostic Catheter	Microcatheter	Diagnostic Catheter
Establishing the Base	 0.018" Detachable Framing Coil	 0.035" Detachable Framing Coil	 0.018" CX Detachable	 0.035" CX Detachable
Filling the Space	 0.018" Detachable HydroCoil <sup>®</sup>	 0.035" Detachable HydroCoil <sup>®</sup>	 0.018" CX Detachable  0.018" Pushable HydroCoil <sup>®</sup>  0.018" Detachable HydroCoil <sup>®</sup>	 0.035" CX Detachable  0.035" Pushable HydroCoil <sup>®</sup>  0.035" Detachable HydroCoil <sup>®</sup>
Minimizing Reperfusion	TERUMO Patented Hydrogel Technology <sup>4</sup>			

\* For every vessel occlusion or aneurysm where a microcatheter or standard diagnostic catheter is required or preferred.

# AZUR<sup>®</sup> CX

Peripheral Coil System

Available in 2 mm and  
3 mm sizes for small  
vessel embolization

## AZUR<sup>®</sup> CX ORDERING INFORMATION

### Detachable 0.018" System/Pack of 1

PRODUCT CODE	LOOP DIAMETER (mm)	LENGTH* (cm)	PRODUCT CODE	LOOP DIAMETER (mm)	LENGTH* (cm)
<b>Expanded portfolio with 2 and 3 mm sizes</b>			45-780828	8	28
45-780202	2	2	45-780928	9	28
45-780204	2	4	45-781032	10	32
45-780304	3	4	45-781238	12	38
45-780308	3	8	45-781434	14	34
45-780413	4	13	45-781639	16	39
45-780516	5	16	45-781836	18	36
45-780620	6	20	45-782040	20	40
45-780724	7	24	—	—	—

### Detachable 0.035" System/Pack of 1

PRODUCT CODE	LOOP DIAMETER (mm)	LENGTH* (cm)	PRODUCT CODE	LOOP DIAMETER (mm)	LENGTH* (cm)
45-750407	4	7	45-750824	8	24
45-750511	5	11	45-751019	10	19
45-750609	6	9	45-751324	13	24
45-750617	6	17	45-751632	16	32
45-750812	8	12	45-752039	20	39

### Catheter ID Requirements

COIL TYPE	INCHES (min-max)	MILLIMETERS (min-max)	REPOSITIONING TIME
AZUR <sup>®</sup> CX 0.018	.019 - .027	.48 - .69	30 minutes
AZUR <sup>®</sup> CX 0.035	.041 - .047	1.04 - 1.19	20 minutes

### Detachable Controller for Use with Detachable Systems/Pack of 5

PRODUCT CODE	PRODUCT DESCRIPTION
45-4001	AZUR <sup>®</sup> Detachment Controller

\* Length is calculated from the tip of the coil to the point of connection with the pusher wire when the coil is straight.

**FIND OUT MORE** ☎ Phone: 800.862.4143 🌐 [terumo.com](http://terumo.com) 📠 Fax: 800.411.5870

#### References:

1. Ding YH, Dai D, Lewis DA, Cloft HJ, Kallmes DF. Angiographic and histologic analysis of experimental aneurysms embolized with platinum coils, Matrix, and HydroCoil. *AJNR Am J Neuroradiol.* 2005 Aug;26(7):1757-63. (animal study)
2. Fanning NF, Berentei Z, Brennan PR, Thornton J. HydroCoil as an adjuvant to bare platinum coil treatment of 100 cerebral aneurysms. *Neuroradiology.* 2006. doi: 10.1007/s00234-006-0166-0. (in-vivo study)
3. Pelage JP. Angiographic and Pathologic Comparison of HydroCoils vs. Fibered Coils Mechanisms of Occlusion and Mid-Term Recanalization in an Animal Model. *GEST.* 2012 (animal study)
4. Plenk H, Killer M, Richling B. Pathophysiologic considerations on HydroCoil- and platinum coil-occluded retrieved human cerebral aneurysms. Presented at ASITN MicroVention Symposium. 2005 (in-vivo study)
5. Data on File. Terumo Corporation. IMS Data.

For Rx only.

Before using refer to **Instructions for Use** for indications, contraindications as well as warnings and precautions at [www.terumo.com](http://www.terumo.com)

©2017 Terumo Medical Corporation. All rights reserved. All brand names are trademarks or registered trademarks of their respective owners. TIS-599-08182017/GGG

**TERUMO**  
INTERVENTIONAL  
SYSTEMS